

FHL3 Antibody (R33604)

Catalog No.	Formulation	Size
R33604-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2275
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This FHL3 antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with FHL3 antibody at 1ug/ml.
Predicted molecular weight: ~31 kDa, can also be observed at 33-34 kDa.

Description

Additional name(s) for this target protein: Four and a half LIM domains 3

Application Notes

Optimal dilution of the FHL3 antibody should be determined by the researcher.

Immunogen

Amino acids DQVLCQGCSQAGP were used as the immunogen for this FHL3 antibody.

Storage

Aliquot and store the FHL3 antibody at -20oC.