


## Flotillin 1 Antibody (Biotin Conjugate) (R35932BTN)

Catalog No.	Formulation	Size
R35932BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10211
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This Flotillin 1 antibody is available for research use only.



Western blot testing of K562 lysate with biotinylated Flotillin 1 antibody at 1ug/ml.  
Expected/observed molecular weight ~49 kDa.

### Description

Additional name(s) for this target protein: FLOT1

### Application Notes

Optimal dilution of the Flotillin 1 antibody should be determined by the researcher.

## Immunogen

Amino acids SISQVNHKPLRTA were used as the immunogen for this Flotillin 1 antibody.

## Storage

Aliquot and store the Flotillin 1 antibody at -20oC.