


FOXG1 Antibody (R34258)

| Catalog No. | Formulation | Size |
|--------------|--|--------|
| R34258-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$ | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|--------------------|---|
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 2290 |
| Applications | Western Blot: 0.3-1ug/ml ELISA (peptide) LOD: 1:64000 |
| Limitations | This FOXG1 antibody is available for research use only. |



Description

Additional name(s) for this target protein: Forkhead box G1, BF1, BF2, FHKL3, FKH2

Application Notes

Optimal dilution of the FOXG1 antibody should be determined by the researcher.

Immunogen

Amino acids HQNQGSSSNPLIH were used as the immunogen for this FOXG1 antibody.

Storage Aliquot and store the FOXG1 antibody at -20oC.