


FOXN1 Antibody (R36106)

Catalog No.	Formulation	Size
R36106-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8456
Applications	Immunohistochemistry (FFPE) : 20-30ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FOXN1 antibody is available for research use only.



IHC testing of FFPE human thymus medulla with FOXN1 antibody at 20ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Forkhead box protein N1

Application Notes

Optimal dilution of the FOXN1 antibody should be determined by the researcher.

Immunogen

Amino acids SVYLSPSSKPVALA were used as the immunogen for this FOXN1 antibody.

Storage Aliquot and store the FOXN1 antibody at -20oC.