


## GAB2 Antibody (R35496)

Catalog No.	Formulation	Size
R35496-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9846
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This GAB2 antibody is available for research use only.



Western blot testing of K562 cell lysate with GAB2 antibody at 0.5ug/ml. Predicted molecular weight: ~75kDa.

## Description

Additional name(s) for this target protein: GRB2-associated binding protein 2

## Application Notes

Optimal dilution of the GAB2 antibody should be determined by the researcher.

## Immunogen

Amino acids QVDKEKTQALQN were used as the immunogen for this GAB2 antibody.

## **Storage**

Aliquot and store the GAB2 antibody at -20oC.