


GATA5 Antibody (R33856)

Catalog No.	Formulation	Size
R33856-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	140628
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GATA5 antibody is available for research use only.



Western blot testing of HepG2 cell lysate with GATA5 antibody at 0.3ug/ml. Predicted molecular weight ~41kDa.

Description

Additional name(s) for this target protein: GATA binding protein 5

Application Notes

Optimal dilution of the GATA5 antibody should be determined by the researcher.

Immunogen

Amino acids ELAARPGWAQTATAD were used as the immunogen for this GATA5 antibody.

Storage

Aliquot and store the GATA5 antibody at -20oC.