


## GATA6 Antibody (F49005)

Catalog No.	Formulation	Size
F49005-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F49005-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

### Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q92908
Applications	Western Blot : 1:1000
Limitations	This GATA6 antibody is available for research use only.



## Description

GATA6 is thought to be important for regulating terminal differentiation and/or proliferation.

## Application Notes

Titration of the GATA6 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 254-279 from the human protein was used as the immunogen for this GATA6 antibody.

## **Storage**

Aliquot the GATA6 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.