


GATAD2B Antibody (F43907)

Catalog No.	Formulation	Size
F43907-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43907-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8WXI9
Applications	Western Blot : 1:1000
Limitations	This GATAD2B antibody is available for research use only.



Description

GATAD2B has transcriptional repressor activity. Targets MBD3 to discrete loci in the nucleus.

Application Notes

Titration of the GATAD2B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1-30 from the human protein was used as the immunogen for this GATAD2B antibody.

Storage

Aliquot the GATAD2B antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.