


GNL3 Antibody (R33842)

Catalog No.	Formulation	Size
R33842-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q9BVP2
Gene ID	26354
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This GNL3 antibody is available for research use only.



Western blot testing of A) HeLa, B) HepG2, and C) K562 nuclear lysates with GNL3 antibody at 2ug/ml. Predicted molecular weight: ~62/60kDa (isoforms 1/2). The observed ~75kDa band likely represents the acetylated form of the protein.

Description

Additional name(s) for this target protein: Guanine nucleotide binding protein-like 3

Application Notes

Optimal dilution of the GNL3 antibody should be determined by the researcher.

Immunogen

Amino acids EKNNAQSIRAIKGPH were used as the immunogen for this GNL3 antibody. **Storage** Aliquot and store the GNL3 antibody at -20oC.