


GRAMD3 Antibody (R35490)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35490-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse |
| Predicted Reactivity | Cow, Dog, Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 65983 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000 |
| Limitations | This GRAMD3 antibody is available for research use only. |



Western blot testing of human temporal cortex lysate with GRAMD3 antibody at 1ug/ml.
Predicted molecular weight: ~48kDa.

Description

Additional name(s) for this target protein: GRAM domain containing 3

Application Notes

Optimal dilution of the GRAMD3 antibody should be determined by the researcher.

Immunogen

Amino acids QTPESSENSRD were used as the immunogen for this GRAMD3 antibody.

Storage

Aliquot and store the GRAMD3 antibody at -20oC.