


## GRIK3 Antibody (R34481)

Catalog No.	Formulation	Size
R34481-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2899
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This GRIK3 antibody is available for research use only.



Western blot testing of rat brain lysate with GRIK3 antibody at 0.3ug/ml. The expected ~110kDa band and additional ~39kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: Glutamate receptor, ionotropic, kainate 3; GLUR7

## Application Notes

Optimal dilution of the GRIK3 antibody should be determined by the researcher.

## Immunogen

Amino acids KIRQLPIDSDDSRP were used as the immunogen for this GRIK3 antibody.

## **Storage**

Aliquot and store the GRIK3 antibody at -20oC.