


## Gsta5 Antibody (R33881)

Catalog No.	Formulation	Size
R33881-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	14859
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Gsta5 antibody is available for research use only.



Western blot testing of rat liver lysate with Gsta5 antibody at 0.3ug/ml. Predicted molecular weight: ~25kDa.

## Description

Additional name(s) for this target protein: Glutathione S-transferase alpha-5

## Application Notes

Optimal dilution of the Gsta5 antibody should be determined by the researcher.

1. This Gsta5 antibody will also detect the homologous mouse Gsta3.

## Immunogen

Amino acids EMDPGIVDNFLLK were used as the immunogen for this Gsta5 antibody.

## Storage

Aliquot and store the Gsta5 antibody at -20oC.