


# **HES6 Antibody (F51760)**

Catalog No.	Formulation	Size
F51760-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F51760-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q96HZ4
Applications	Western Blot: 1:1000
Limitations	This HES6 antibody is available for research use only.



#### **Description**

HES6 is a member of a subfamily of basic helix-loop-helix transcription repressors that have homology to the Drosophila enhancer of split genes. The protein functions as a cofactor, interacting with other transcription factors through a tetrapeptide domain in its C-terminus.

## **Application Notes**

Titration of the HES6 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### **Immunogen**

A portion of amino acids 6-34 from the human protein was used as the immunogen for this HES6 antibody. **Storage** Aliquot the HES6 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com Copyright  $\ensuremath{\texttt{@}}$  NSJ Bioreagents. All rights reserved.