


# **HLA-DQA2** Antibody (R35277)

Catalog No.	Formulation	Size
	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3118; 3117
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HLA-DQA2 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Major histocompatibility complex, class II, DQ alpha 2

#### **Application Notes**

Optimal dilution of the HLA-DQA2 antibody should be determined by the researcher.

1. This HLA-DQA2 antibody will also detect the highly homologous HLA-DQA1.

## **Immunogen**

Amino acids QGLRSVGASRH were used as the immunogen for this HLA-DQA2 antibody.(1)

# **Storage**

Aliquot and store the HLA-DQA2 antibody at -20oC.