


HLA-DQA2 Antibody (R35457)

Catalog No.	Formulation	Size
R35457-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3118
Applications	Western Blot: 0.1-0.3 ug/ml ELISA (peptide) LOD: 1:16000
Limitations	This HLA-DQA2 antibody is available for research use only.



Description

Additional name(s) for this target protein: Major histocompatibility complex, class II, DQ alpha 2

Application Notes

Optimal dilution of the HLA-DQA2 antibody should be determined by the researcher.

Immunogen

Amino acids QSHGPSGQYTH were used as the immunogen for this HLA-DQA2 antibody.

Storage Aliquot and store the HLA-DQA2 antibody at -20oC.