


Hras1 Antibody (F47429)

Catalog No.	Formulation	Size
F47429-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47429-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Human, Rat, Chicken
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q61411
Applications	Western Blot : 1:1000
Limitations	This Hras1 antibody is available for research use only.



Hras1 antibody western blot analysis in mouse bladder tissue lysate. Predicted molecular weight ~ 21 kDa.

Description

Ras proteins bind GDP/GTP and possess intrinsic GTPase activity.

Application Notes

Titration of the Hras1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 151-179 from the mouse protein was used as the immunogen for this Hras1 antibody.

Storage

Aliquot the Hras1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.