


HSD11B1 Antibody (R33358)

Catalog No.	Formulation	Size
R33358-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3290
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This HSD11B1 antibody is available for research use only.



Western blot testing of human liver lysate with HSD11B1 antibody at 0.3ug/ml. Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: 11-beta-hydroxysteroid dehydrogenase 1

Application Notes

Optimal dilution of the HSD11B1 antibody should be determined by the researcher.

Immunogen

Amino acids CTSYNMDRFINK were used as the immunogen for this HSD11B1 antibody.

Storage

Aliquot and store the HSD11B1 antibody at -20oC.