


Hsp31p Antibody (R34015)

Catalog No.	Formulation	Size
R34015-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Citations (1)	Bulk quote request
Citations (1)	Bulk quote reques

Availability	1-3 business days
Species Reactivity	S. cerevisiae
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q04432
Applications	Western Blot: 0.05-0.2ug/ml ELISA (peptide) LOD: 1:128000
Limitations	This Hsp31p antibody is available for research use only.



Western blot testing of yeast lysate (wt) and the KO (delta) in the left panel and the Ponceau stain in the right panel. Hsp31p antibody was used at 0.05ug/ml.

Description

Additional name(s) for this target protein: Heat shock protein 31

Application Notes

Optimal dilution of the Hsp31p antibody should be determined by the researcher.

Immunogen

Amino acids ETGKFGWDEHSLAKD were used as the immunogen for this Hsp31p antibody.

Storage Aliquot and store the Hsp31p antibody at -20oC.						