


IA2 Antibody (R35270)

Catalog No.	Formulation	Size
R35270-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5798
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This IA2 antibody is available for research use only.



Description

Additional name(s) for this target protein: Islet cell antigen 512; PTPRN; protein tyrosine phosphatase receptor type N

Application Notes

Optimal dilution of the IA2 antibody should be determined by the researcher.

Immunogen

Amino acids QGEGNIKKNRHPD were used as the immunogen for this IA2 antibody. **Storage** Aliquot and store the IA2 antibody at -20oC.