


# ICAM4 Antibody (R35626)

Catalog No.	Formulation	Size
R35626-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3386
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ICAM4 antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: Intercellular adhesion molecule 4

#### **Application Notes**

Optimal dilution of the ICAM4 antibody should be determined by the researcher.

1. This ICAM4 antibody will detect isoforms 1 and 2.

# Immunogen

Amino acids RHGSRVIYSESLER were used as the immunogen for this ICAM4 antibody.

### **Storage**

Aliquot and store the ICAM4 antibody at -20oC.