


IL13RA2 Antibody (R33538)

Catalog No.	Formulation	Size
R33538-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3598
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IL13RA2 antibody is available for research use only.



Western blot testing of HeLa lysate with IL13RA2 antibody at 1ug/ml. Predicted molecular weight: ~44/60kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Interleukin 13 receptor alpha 2

Application Notes

Optimal dilution of the IL13RA2 antibody should be determined by the researcher.

Immunogen

Amino acids DLNKGIEAKIHT were used as the immunogen for this IL13RA2 antibody.

Storage

Aliquot and store the IL13RA2 antibody at -20oC.

References (1)