


IL17RD Antibody (F51826)

Catalog No.	Formulation	Size
F51826-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F51826-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8NFM7
Applications	Western Blot : 1:1000
Limitations	This IL17RD antibody is available for research use only.



Western blot analysis of IL17RD antibody and MDA-MB231 lysate.

Description

Fibroblast growth factors are secreted proteins involved in cellular proliferation, migration, differentiation, and survival. FGF activity is negatively regulated by members of the 'sprouty' family. SEF/IL17RD is a modulator of FGF signaling.

Application Notes

Titration of the IL17RD antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 217-245 from the human protein was used as the immunogen for this IL17RD antibody.

Storage

Aliquot the IL17RD antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.