


IL17RE Antibody (F47506)

Catalog No.	Formulation	Size
F47506-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47506-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8NFR9
Applications	Western Blot : 1:1000
Limitations	This IL17RE antibody is available for research use only.



Description

This gene encodes a protein that shares limited similarity with the interleukin-17 receptor. The function of the encoded protein has not yet been determined. Alternatively spliced transcript variants have been described.

Application Notes

Titration of the IL17RE antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 87-116 from the human protein was used as the immunogen for this IL17RE antibody.

Storage

Aliquot the IL17RE antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.