


Irak4 Antibody (F43563)

Catalog No.	Formulation	Size
F43563-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43563-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8R4K2
Applications	Western Blot: 1:1000
Limitations	This Irak4 antibody is available for research use only.



Description

Irak4 is required for the efficient recruitment of IRAK1 to the IL-1 receptor complex following IL-1 engagement, triggering intracellular signaling cascades leading to transcriptional up-regulation and mRNA stabilization. Phosphorylates IRAK1 (By similarity).

Application Notes

Titration of the Irak4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 422-449 from the mouse protein was used as the immunogen for this Irak4 antibody. **Storage** Aliquot the Irak4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com Copyright $\ensuremath{\texttt{@}}$ NSJ Bioreagents. All rights reserved.