


## Itih4 Antibody (R34563)

Catalog No.	Formulation	Size
R34563-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	16427
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Itih4 antibody is available for research use only.



Western blot testing of rat liver lysate with Itih4 antibody at 0.03ug/ml. Predicted molecular weight: ~103/150kDa (unmodified/glycosylated).

## Description

Additional name(s) for this target protein: Inter-alpha (globulin) inhibitor H4

## Application Notes

Optimal dilution of the Itih4 antibody should be determined by the researcher.

## Immunogen

Amino acids KPEGQEQFQVAEKP were used as the immunogen for this Itih4 antibody.

## Storage

Aliquot and store the Itih4 antibody at -20oC.