


KIAA0319 Antibody (R35274)

Catalog No.	Formulation	Size
R35274-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	210108
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This KIAA0319 antibody is available for research use only.



Description

Additional name(s) for this target protein: D130043K22Rik, 4930451E12Rik

Application Notes

Optimal dilution of the KIAA0319 antibody should be determined by the researcher.

Immunogen

Amino acids QGKIKQENKPTLH were used as the immunogen for this mouse KIAA0319 antibody. **Storage** Aliquot and store the KIAA0319 antibody at -20oC.