


KPNA3 Antibody (Biotin Conjugate) (R33584BTN)

Catalog No.	Formulation	Size
R33584BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3839
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This KPNA3 antibody is available for research use only.



Western blot testing of rat testis lysate with biotinylated KPNA3 antibody at 1 ug/ml.
Predicted molecular weight ~58kDa.

Description

Additional name(s) for this target protein: Karyopherin alpha 3; IPOA4; importin alpha 4

Application Notes

Optimal dilution of the KPNA3 antibody should be determined by the researcher.

Immunogen

Amino acids DPTANLQTKEFNF were used as the immunogen for this KPNA3 antibody.

Storage

Aliquot and store the KPNA3 antibody at -20oC.