


Lingo1 Antibody (F54318)


Catalog No.	Formulation	Size
F54318-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54318-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity purified
UniProt	Q9D1T0
Applications	Western Blot : 1:500-1:2000 Immunohistochemistry (FFPE) : 1:25
Limitations	This Lingo1 antibody is available for research use only.



IHC testing of FFPE human brain tissue with Lingo1 antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.



IHC testing of FFPE human brain tissue with Lingo1 antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.



Description

'Leucine-rich repeat and immunoglobulin domain-containing protein 1' is the functional component of the Nogo receptor signaling complex (RTN4R/NGFR) in RhoA activation responsible for some inhibition of axonal regeneration by myelin-associated factors. LINGO1 is also an important negative regulator of oligodentrocyte differentiation and axonal myelination. [UniProt]

Application Notes

The stated application concentrations are suggested starting points. Titration of the Lingo1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 575-603 from the mouse protein was used as the immunogen for the Lingo1 Antibody.

Storage

Aliquot the Lingo1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.