


LOC147646 Antibody (F40018)

Catalog No.	Formulation	Size
F40018-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40018-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	G1R1R7
Applications	Western Blot : 1:1000
Limitations	This LOC147646 antibody is available for research use only.



Description

The function of this protein remains unknown.

Application Notes

Titration of the LOC147646 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 147-176 from the nomle protein was used as the immunogen for this LOC147646 antibody.

Storage

Aliquot the LOC147646 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.