


# MAFA Antibody (F44543)

Catalog No.	Formulation	Size
F44543-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F44543-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8NHW3
Applications	Western Blot : 1:1000
Limitations	This MAFA antibody is available for research use only.



# **Description**

MAFA is a transcription factor that binds RIPE3b, a conserved enhancer element that regulates pancreatic beta cell-specific expression of the insulin gene.

## **Application Notes**

Titration of the MAFA antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### **Immunogen**

A portion of amino acids 286-315 from the human protein was used as the immunogen for this MAFA antibody.

# **Storage**

Aliquot the MAFA antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.