


MAFA Antibody (F51722)

Catalog No.	Formulation	Size
F51722-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F51722-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q8NHW3
Applications	Western Blot : 1:1000
Limitations	This MAFA antibody is available for research use only.



Western blot analysis of MAFA antibody and mouse spleen tissue lysate. Predicted molecular weight ~37 kDa.

Description

MAFA is a transcription factor that binds RIPE3b, a conserved enhancer element that regulates pancreatic beta cell-specific expression of the insulin gene.

Application Notes

Titration of the MAFA antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 203-231 from the human protein was used as the immunogen for this MAFA antibody.

Storage

Aliquot the MAFA antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.