


MCT-1 Antibody (F49765)

Catalog No.	Formulation	Size
F49765-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F49765-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Bovine, Chicken, Xenopus, Zebrafish
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q9ULC4
Applications	Western Blot : 1:1000
Limitations	This MCT-1 antibody is available for research use only.



Description

MCTS1 play a role in cell cycle regulation; decreases cell doubling time and anchorage-dependent growth; shortens the duration of G1 transit time and G1/S transition.

Application Notes

Titration of the MCT-1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1-30 from the human protein was used as the immunogen for this MCT-1 antibody.

Storage

Aliquot the MCT-1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.