


MPP5 Antibody (R34832)

Catalog No.	Formulation	Size
R34832-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64398
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MPP5 antibody is available for research use only.



Description

Additional name(s) for this target protein: Membrane protein palmitoylated 5, PALS, PALS1; MPP5

Application Notes

Optimal dilution of the MPP5 antibody should be determined by the researcher.

Immunogen

Amino acids NKASPPFPLISNAQ were used as the immunogen for this MPP5 antibody. **Storage** Aliquot and store the MPP5 antibody at -20oC.