


MSH5 Antibody (R34074)

Catalog No.	Formulation	Size
R34074-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4439
Applications	Immunohistochemistry (FFPE) : 2.5-5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This MSH5 antibody is available for research use only.



IHC testing of FFPE human testis with MSH5 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: MutS homolog 5

Application Notes

Optimal dilution of the MSH5 antibody should be determined by the researcher.

Immunogen

Amino acids GDLHCEIRDQET were used as the immunogen for this MSH5 antibody. **Storage** Aliquot and store the MSH5 antibody at -20oC.