


# NIrp12 Antibody (R33401)

Catalog No.	Formulation	Size
R33401-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	378425
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Nlrp12 antibody is available for research use only.



# **Description**

Additional name(s) for this target protein: NLR family, pyrin domain containing 12

#### **Application Notes**

Optimal dilution of the Nlrp12 antibody should be determined by the researcher.

## Immunogen

Amino acids DDPPEPSGVQTQST were used as the immunogen for this Nlrp12 antibody.

Storage Aliquot and store the NIrp12 antibody at -20oC.						