


# NNMT Antibody (R34976)

Catalog No.	Formulation	Size
R34976-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4837
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This NNMT antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Nicotinamide N-methyltransferase

### **Application Notes**

Optimal dilution of the NNMT antibody should be determined by the researcher.

### Immunogen

Amino acids PDLPTYCRALRN were used as the immunogen for this NNMT antibody.

# **Storage** Aliquot and store the NNMT antibody at -20oC.