


# NRG3 Antibody (R35689)

Catalog No.	Formulation	Size
R35689-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10718
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This NRG3 antibody is available for research use only.



#### **Description**

Additional name(s) for this target protein: Neuregulin 3

## **Application Notes**

Optimal dilution of the NRG3 antibody should be determined by the researcher.

### **Immunogen**

Amino acids RNEIQRDSALTK were used as the immunogen for this NRG3 antibody. **Storage** Aliquot and store the NRG3 antibody at -20oC.