


# OAS2 Antibody (R35659)

Catalog No.	Formulation	Size
R35659-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4939
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This OAS2 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: 2'-5'-oligoadenylate synthase 2

## **Application Notes**

Optimal dilution of the OAS2 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids RKTVLRGNSDGT were used as the immunogen for this OAS2 antibody. This sequence is common to

isoforms 1 (82 kDa, can be observed at 71kDa) and 3 (~20kDa).

Storage

Aliquot and store the OAS2 antibody at -20oC.