


Olig3 Antibody (F53071)

Catalog No.	Formulation	Size
F53071-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53071-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q7RTU3
Applications	Western Blot : 1:1000
Limitations	This Olig3 antibody is available for research use only.



Western blot analysis of lysate from (1) mouse brain and (2) mouse skeletal muscle tissue lysate using Olig3 antibody at 1:1000.

Description

May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons.

Application Notes

Titration of the Olig3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This Olig3 antibody was produced from a rabbit immunized with a KLH conjugated synthetic peptide between 15-48

amino acids from the N-terminal region of human OLIG3.

Storage

Aliquot the Olig3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.