


# OR1B1 Antibody (F54215)

Catalog No.	Formulation	Size
F54215-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54215-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal
Isotype	Rabbit Ig
Purity	Peptide affinity purified
UniProt	Q8NGR6
Applications	Western Blot: 1:500
Limitations	This OR1B1 antibody is available for research use only.



## **Description**

Olfactory receptor 1B1 is an odorant receptor.

# **Application Notes**

The stated application concentrations are suggested starting points. Titration of the OR1B1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### **Immunogen**

A portion of amino acids 149-183 from the human protein was used as the immunogen for the OR1B1 antibody.

# **Storage** Aliquot the OR1B1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.