


P130Cas Antibody (R34605)


Catalog No.	Formulation	Size
R34605-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9564
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This P130Cas antibody is available for research use only.



Western blot testing of MCF7 [A] and human breast cancer [B] lysates with P130Cas antibody at 1ug/ml. Predicted molecular weight: ~98kDa.



IHC staining of FFPE human uterus with P130Cas antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Breast cancer estrogen resistance 1; BCAR1

Application Notes

Optimal dilution of the P130Cas antibody should be determined by the researcher.

1. This P130Cas antibody will detect isoforms 1/2/3/5/6/7.

Immunogen

Amino acids KQTSTFSKQTPHH were used as the immunogen for this P130Cas antibody.

Storage

Aliquot and store the P130Cas antibody at -20oC.