


PGA3 Antibody (F40034)

Catalog No.	Formulation	Size
F40034-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40034-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P0DJD8
Applications	Western Blot: 1:1000
Limitations	This PGA3 antibody is available for research use only.



Description

PGA3 shows particularly broad specificity; although bonds involving phenylalanine and leucine are preferred, many others are also cleaved to some extent.

Application Notes

Titration of the PGA3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 45-73 from the human protein was used as the immunogen for this PGA3 antibody.

Storage Aliquot the PGA3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.