

PGA4 Antibody (F40025)


Catalog No.	Formulation	Size
F40025-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40025-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Rabbit
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P0DJD7
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This PGA4 antibody is available for research use only.



PGA4 antibody western blot analysis in WiDr lysate.



PGA4 antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human stomach tissue.

Description

PGA4 shows particularly broad specificity; although bonds involving phenylalanine and leucine are preferred, many others are also cleaved to some extent.

Application Notes

Titration of the PGA4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 353-381 from the human protein was used as the immunogen for this PGA4 antibody.

Storage

Aliquot the PGA4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.