


PGA5 Antibody (F40041)

Catalog No.	Formulation	Size
F40041-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40041-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Rabbit
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P0DJD9
Applications	Western Blot : 1:1000
Limitations	This PGA5 antibody is available for research use only.



Description

PGA3 shows particularly broad specificity; although bonds involving phenylalanine and leucine are preferred, many others are also cleaved to some extent.

Application Notes

Titration of the PGA5 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 354-383 from the human protein was used as the immunogen for this PGA5 antibody.

Storage

Aliquot the PGA5 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.