


Phospho-Bad Antibody (pS97) (F48679)

Catalog No.	Formulation	Size
F48679-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F48679-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	O35147
Applications	Dot Blot : 1:500
Limitations	This phospho-Bad antibody is available for research use only.



Description

Bad may be involved in induction of programmed cell death [RGD].

Application Notes

Titration of the phospho-Bad antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This phospho-Bad antibody was produced from rabbits immunized with a KLH conjugated synthetic phosphopeptide

corresponding to amino acid residues surrounding pS97 of rat Bad. **Storage** Aliquot the phospho-Bad antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com Copyright $\ensuremath{\texttt{@}}$ NSJ Bioreagents. All rights reserved.