


## Pleckstrin Antibody (R34045)

Catalog No.	Formulation	Size
R34045-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5341
Applications	Western Blot : 0.05-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Pleckstrin antibody is available for research use only.



Western blot testing of human PBMC lysate with Pleckstrin antibody at 0.05ug/ml.  
Predicted molecular weight: ~40kDa, observed here at ~45kDa.

## Description

Additional name(s) for this target protein: P47; PLEK

## Application Notes

Optimal dilution of the Pleckstrin antibody should be determined by the researcher.

## Immunogen

Amino acids EWIKAIQMASRTGK were used as the immunogen for this Pleckstrin antibody.

## Storage

Aliquot and store the Pleckstrin antibody at -20oC.