


Prdm9 Antibody (R35239)

Catalog No.	Formulation	Size
R35239-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	365155
Applications	Western Blot : 2-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Prdm9 antibody is available for research use only.



Description

Additional name(s) for this target protein: PR domain containing 9

Application Notes

Optimal dilution of the Prdm9 antibody should be determined by the researcher.

Immunogen

Amino acids QEHFDSQNKNDK were used as the immunogen for this Prdm9 antibody.

Storage Aliquot and store the Prdm9 antibody at -20oC.