


PRMT5 Antibody (R35528)


Catalog No.	Formulation	Size
R35528-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10419
Applications	Immunofluorescence : 10ug/ml Flow Cytometry : 10ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PRMT5 antibody is available for research use only.



IF/ICC staining of fixed and permeabilized human U-2 OS cells with PRMT5 antibody (green) at 10ug/ml and DAPI nuclear stain (blue).



IF/ICC staining of fixed and permeabilized human A431 cells with PRMT5 antibody (green) at 10ug/ml and DAPI nuclear stain (blue).



FACS testing of fixed and permeabilized human A431 cells with PRMT5 antibody (blue) at 10ug/ml and naive goat lg (black).

Description

Additional name(s) for this target protein: Protein arginine methyltransferase 5

Application Notes

Optimal dilution of the PRMT5 antibody should be determined by the researcher.

Immunogen

Amino acids RACREKDRDPEAQ were used as the immunogen for this PRMT5 antibody.

Storage

Aliquot and store the PRMT5 antibody at -20oC.