


Rab17 Antibody (R35645)

Catalog No.	Formulation	Size
R35645-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	19329
Applications	Western Blot : 0.2-0.6ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Rab17 antibody is available for research use only.



Description

Additional name(s) for this target protein: Ras related protein Rab17

Application Notes

Optimal dilution of the Rab17 antibody should be determined by the researcher.

Immunogen

Amino acids RKDSFHKAQQWLED were used as the immunogen for this Rab17 antibody. **Storage** Aliquot and store the Rab17 antibody at -20oC.