


# RAB7B Antibody (F47662)

Catalog No.	Formulation	Size
F47662-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47662-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q96AH8
Applications	Western Blot: 1:1000
Limitations	This RAB7B antibody is available for research use only.



# **Description**

Involved in protein transport. Likely involved in vesicular traffic (By similarity). [UniProt]

# **Application Notes**

Titration of the RAB7B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

# **Immunogen**

A portion of amino acids 168-197 from the human protein was used as the immunogen for this RAB7B antibody.

# **Storage** Aliquot the RAB7B antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.